

3D Sugical Mesh

- Used for treatment of inguinal hernia repairment.

Soft material


Less foreign body sensation after implantation and high comfort

Low shrinkage

Less pain after implantation

Three-dimensional stereotyping

Easy positioning, anti-displacement, anti-recurrence



Enhanced sealed Edge Design

Easy to unfold, tear-proof, wrinkle-proof


Large aperture

Special woven structure, which is conducive to the rapid growth of tissue

Recommedation Spec

Shape	Spec	Width (±8mm)	Length (±10mm)	Height (±6mm)
Left mesh	BSLB79×134	79	134	18
Left mesh	BSLB103×157	103	157	21.5
Right mesh	BSRB79×134	79	134	18
Right mesh	BSRB103×157	103	157	21.5

3D Sugical Mesh



Lightly weight suitable for hernia
Perfectly shape easier to place



Sinolinks Medical Innovation, Inc.
Add.: No.10 Ziwei Road, Zhonglou Zone, Changzhou Jiangsu, China, 213023
Tel: 0519-88021690
E-mail: sales@sinolinks.com.cn
Web: www.sinolinks.com.cn