



Excellent sealing, smooth insertion

Unique Four-valve plus ultra-thin sealing membrane design, to ensure tightness and reduce the frictional resistance of instruments passing through.



Detachable Housing

Simply press the knob to eject the seal assembly. Quick and easy to operate with one hand.



Better visualization

Visible light-colored puncture rod for better light transmission, convenient to see the entire process of puncturing abdominal wall. Avoiding damage to the abdominal tissue.

Ordering information

Code	Diameter	working length
AVTC 5 - 100	5	100
AVTC 10 - 100	10	100
AVTC 12 - 100	12	100
AVTC 15 - 100	15	100

- Sinolinks Medical Innovation, Inc.
- Add.: No.10 Ziwei Road, Zhonglou Zone, Changzhou, Jiangsu, China, 213023
- Tel:0519-88021690
- E-mail:sales@sinolinks.com.cn
- Web site:www.sinolinks.com.cn

Disposable Trocar



Application

Used for minimally invasive surgical procedures to establish a path of entry for endoscopic instruments for primary and secondary insertions.